

Affymetrix Genome-Wide SNP Array 6.0

Affymetrix Genome-Wide SNP Array 6.0 产品的特点是超过 1,800,000 个遗传变异标志物,包括超过 906,600 个 SNP 和超过 946,000 个用于检测拷贝数变化的探针。大约有 482,000 个 SNPs 来自于前代产品 500K 和 SNP5.0 芯片。剩下 424,000s 个 SNP 包括了来源于国际 HapMap 计划中的标签 SNPs, X, Y 染色体和线粒体上更具代表性的 SNPs,以及来自于重组热点区域和 500K 芯片设计完成后新加入 dbSNP 数据库的 SNPs。SNP6.0 还包括了 202,000 个用于检测 5,677 个已知拷贝数变异区域的探针,这些区域来源的于多伦多基因组变异体数据库。该数据库中的每个 3,182 个非重叠片断区域分别平均用 61 个探针来检测。除了检测这些已知的拷贝数多态区域,还有超过 744,000 个探针平均分配到整个基因组上,用来发现求知的拷贝数变异区域。设计方面 SNP 6.0 芯片只取最佳 SNP 位点和每组等到位基因的一个完全匹配,并有三次重复检测,所以总共只有六个探针。

具体参数如下:

超过 906,600 个 SNP, 包括:

Unbiased selection of 482,000 SNPs; historical SNPs from the SNP Array 5.0

Selection of additional 424,000 SNPs

Tag SNPs

SNPs from chromosomes X and Y

Mitochondrial SNPs

New SNPs added to the dbSNP database

SNPs in recombination hotspots

超过 946,000 的用于检测拷贝数变化的探针, 包括:

202,000 probes targeting 5,677 known CNV regions from the Toronto Database of Genomic Variants

Regions resolve into 3,182 distinct, non-overlapping segments; on average 61 probes per region
744,000 probes,

evenly spaced along the genome

样品要求

1. 样品纯度：OD 260/280 值应在 1.7~2.0 之间；RNA 应该去除干净；不得有其它个体或其它物种的 DNA 污染。
2. 样品浓度：浓度不低于 55ng/μl。
3. 样品总量：每个样品总量不少于 1μg。
4. 样品溶剂：溶解在 Reduced TE (10mM Tris, pH 8.0, 0.1mM EDTA)中。
5. 样品运输：DNA 低温运输 (-20℃)；在运输过程中请用 parafilm 将管口密封好，以防出现污染。

试验流程

1. DNA抽提；
2. 样品质检：用琼脂糖凝胶电泳或LAB-ON-CHIP系统对样品DNA进行质量检测和定量；
3. 实验过程：
 - a、扩增，标记；
 - b、杂交、洗脱（具体见左图）；
4. 图像扫描：用Affymetrix Scanner3000 7G激光共聚焦扫描仪对杂交结束后的芯片进行扫描；
5. 数据处理：确定检出的SNP位点相应的核苷酸。

数据分析流程

1. SNP 位点过滤
 - 1) Nocal rate \geq 10%
 - 2) 最小等位基因频率，Minor allele frequency $<$ 0.05
 - 3) 不符合哈维平衡的位点，HWE P value $<$ 0.001 的位点

2. 频数统计：计算 case/control 组中各等位基因和基因型的频率
3. 关联分析
 - 1) 样品群的分层分析
 - 2) 卡方检验、Fisher 精确检验，比较各 allele 的频率分布在 case 和 control 两组中是否具有统计学差异
 - 3) 卡方检验、Fisher 精确检验和 Cochran-Armitage 趋势检验,比较各基因型频率分布在两组中是否具有统计学差异
4. 单倍型分析
 - 1) 选取包含显著的 SNP 的区段进行单倍型分析，以进一步定位候选区段
 - 2) 频率较高的几种单倍型各自的频率，在 case/control 两组中的分布状况，卡方统计量及 p 值，判断单倍型与疾病相关性

结果报告

- 1、实验数据：原始数据（cel 文件），基因芯片数据总表，SNP 位点信息。
- 2、实验文件：SNP 芯片实验操作方法。
- 3、实验图像：SNP 芯片扫描图（jpg 格式），所有样品的质检电泳图。
- 4、质检文件：DNA 质检报告，SNP Call rate。

Introduction

The new Affymetrix? Genome-Wide Human SNP Array 6.0 contains more than 906,600 single nucleotide polymorphisms (SNPs) and more than 946,000 probes for the detection of copy number variation. SNPs on the array are present on 200 to 1,100 base pairs (bp) Nsp I or Sty I digested fragments in the human genome, and are amplified using the Genome-Wide Human SNP Nsp/Sty Assay Kit 5.0/6.0. This assay, which is also compatible with the SNP Array 5.0, now combines the Nsp and Sty fractions previously assayed on two separate arrays.

SNPs on the Genome-Wide Human SNP Array 6.0 were screened in more than 500 distinct samples, including 270 HapMap samples and separate diversity samples. Approximately 482,000 SNPs are derived from the previous-generation Mapping 500K and SNP 5.0 Arrays. The remaining 424,000 SNPs include tag SNP markers derived from the International HapMap Project, better representation of SNPs on chromosomes X and Y, mitochondrial SNPs, SNPs in recombination hotspots and new SNPs added to the dbSNP database after completion of the Mapping 500K Array. The array also contains 202,000 probes targeting 5,677 known regions of copy number variation from the Toronto Database of Genomic Variants. These regions resolve into 3,182 distinct, nonoverlapping segments, each interrogated with an average of 61 probes. In addition to the interrogation of these regions of known copy number polymorphism, more than 744,000 probes were chosen, evenly spaced along the genome, to enable the detection of novel copy number variation.

The median inter-marker distance taken over all 1.8 million SNP and copy number markers combined is less than 700 bases.

THE WHOLE-GENOME SAMPLING ASSAY

The Affymetrix® Genome-Wide Human SNP Nsp/Sty Assay Kit 5.0/6.0 [P/N 901152, 901015] is validated for use in conjunction with the Genome-Wide Human SNP Array 6.0. Briefly, total genomic DNA (500 ng) is digested with Nsp I and Sty I restriction enzymes and ligated to adaptors that recognize the cohesive 4 bp overhangs. All fragments resulting from restriction enzyme digestion, regardless of size, are substrates for adaptor ligation. A generic primer that recognizes the adaptor sequence is used to amplify adaptor-ligated DNA fragments. PCR conditions have been optimized to preferentially amplify fragments in the 200 to 1,100 bp size range. PCR amplification products for each restriction enzyme digest are combined and purified using polystyrene beads. The amplified DNA is then fragmented, labeled and hybridized to a Genome-Wide Human SNP Array 6.0.

The Genome-Wide Human SNP Nsp/Sty Assay Kit contains validated and qualified reagents for

the most critical steps in the assay. This includes the PCR primer and adaptors, reagents to fragment and label the PCR products and control reagents. Kits are available for either 50 or 100 reactions (see ordering information). Whole-genome-amplified material prepared by the Qiagen REPLI-g? kits may also be used as the starting material for the Genome-Wide Human SNP Assay Kit.

Performance Data

To test the performance of the SNP Array 6.0,


Affymetrix ran the 270 samples from the International HapMap Project. In addition, two external sites ran a plate of 44 HapMap DNAs, which includes 30 unique samples, 10

trios and five samples with multiple replicates. The arrays that passed the QC call rate threshold were analyzed using the Birdseed algorithm at the default setting of 0.1. The average call rate for each set was greater than

99 percent, and the concordance with HapMap genotypes was observed to be greater than or equal to 99.7 percent. For the 10 trios, the Mendelian inheritance consistency was found to be greater than 99.9 percent. Reproducibility was measured at 99.9 percent.

Data analyzed with Birdseed (0.01).			
	270 HapMap	Site 1	Site 2
Call Rate	99.8	99.7	99.7
HapMap Concordance	99.8	99.7	99.8
Mendelian Consistency	99.97	99.95	99.96
Reproducibility	NA	99.9	99.9
SNP Completeness*	99.9	99.7	99.8

*SNP completeness is defined as the proportion of SNPs with per-SNP call rate greater than 85 percent.



GENOTYPE CALLS USING AFFYMETRIX

GENOTYPING CONSOLE

The Genome-Wide Human SNP Array 6.0 is used in conjunction with Affymetrix®

Genotyping Console, which implements a novel genotype calling algorithm called Birdseed.

Birdseed is an evolution of the RLMM genotype calling algorithm¹. It performs a multiple-chip analysis to estimate a signal intensity for each allele of each SNP, fitting probe-specific effects to increase precision (like the BRLMM-P algorithm developed for the Genome-Wide Human SNP Array 5.0). It then makes genotype calls by fitting a Gaussian mixture model in the two-dimensional A-signal vs. B-signal space, using SNP-specific models to improve accuracy.

Genotyping Console was designed to streamline genotyping calling and genotyping quality

control. In addition to the algorithm, features include automated QC that sorts samples by QC call rate, visualization of QC metrics across samples, SNP cluster visualization and signature SNP list for tracking. Seventy-two signature SNPs enable verification of a sample's identity by comparing the genotype calls to different technologies or other known reference. Flexible SNP filter and export are also included to enable downstream analysis. In addition to the Birdseed algorithm, Genotyping Console supports the BRLMMP algorithm for the SNP Array 5.0. Refer to the Affymetrix? Genome-Wide Human SNP Nsp/Sty 5.0 Assay 5.0 or 6.0 Manual (P/N 702419-2, P/N 702054) for details on the QC call rate thresholds, as well as procedures on DNA target preparation, target hybridization, fluidics setup, array scanning and data analysis.

Common variants in DGKK are strongly associated with risk of hypospadias

van der Zanden, L. F. M. *et al.* Nature Genetics 43(1), NPG, 2011 [PubMed](#)

Genome-wide association study identifies susceptibility loci for polycystic ovary syndrome on chromosome 2p16.3, 2p21 and 9q33.3

Chen, Z. J. *et al.* Nature Genetics 43(1), NPG, 2011 [PubMed](#)

Genome-wide association analysis in primary sclerosing cholangitis identifies two non-HLA susceptibility loci

Melum, E. *et al.* Nature Genetics 43(1), 17-9, 2011 [PubMed](#)

Common variants near ATM are associated with glycemic response to metformin in type 2 diabetes

The GoDARTS and UKPDS Diabetes Pharmacogenetics Study Group *et al.* Nature Genetics 43(2), NPG, 2011 [PubMed](#)

Amplification of the 20q Chromosomal Arm Occurs Early in Tumorigenic Transformation and May Initiate Cancer

Tabach, Y. *et al.* Plos One 6(1), NPG, 2011 [PubMed](#)

Genes and biological processes commonly disrupted in rare and heterogeneous developmental delay syndromes

Shaikh, T. H. *et al.* Human Molecular Genetics 20(5), 880-93, 2011 [PubMed](#)

Cohen syndrome diagnosis using whole genome arrays

Rivera-Brugues, N. *et al.* Journal of Medical Genetics 48(2), 136-40, 2011 [PubMed](#)

Tartrate-resistant acid phosphatase deficiency causes a bone dysplasia with autoimmunity and a type I interferon expression signature

Briggs, T. A. *et al.* Nature Genetics 43(2), NPG, 2011 [PubMed](#)

Identification of Novel Susceptibility Loci for Kawasaki Disease in a Han Chinese Population by a Genome-Wide Association Study

Tsai, F. J. *et al.* Plos One 6(2), NPG, 2011 [PubMed](#)

Genome-wide association study identifies a novel susceptibility locus at 6p21.3 among familial CLL

Slager, S. L. *et al.* Blood 117(6), 1911-6, 2011 [PubMed](#)

A New Syndrome With Multiple Capillary Malformations, Intractable Seizures, and Brain and Limb Anomalies

Carter, M. T. *et al.* American Journal of Medical Genetics Part a 155A(2), 301-6, 2011 [PubMed](#)

Intragenic Deletions of IL1RAPL1: Report of Two Cases and Review of the Literature

Behnecke, A. *et al.* American Journal of Medical Genetics Part a 155A(2), 372-9, 2011 [PubMed](#)

47, XY, +der(Y), t(X;Y)(p21.1;p11.2): A Unique Case of XY Sex Reversal

Zarate, Y. A. *et al.* American Journal of Medical Genetics Part a 155A(2), 386-91, 2011 [PubMed](#)

Replication and Functional Genomic Analyses of the Breast Cancer Susceptibility Locus at 6q25.1 Generalize Its Importance in Women of Chinese, Japanese, and European Ancestry

Cai, Q. *et al.* *Cancer Research* 71(4), 1344-55, 2011 [PubMed](#)

High incidence of recurrent copy number variants in patients with isolated and syndromic Mýllerian aplasia

Nik-Zainal, S. *et al.* *Journal of Medical Genetics* 48(3), 197-204, 2011 [PubMed](#)

A Missense Mutation in DHDDS, Encoding Dehydrodolichyl Diphosphate Synthase, Is Associated with Autosomal-Recessive Retinitis Pigmentosa in Ashkenazi Jews

Zelinger, L. *et al.* *American Journal of Human Genetics* 88(2), 207-15, 2011 [PubMed](#)

Bimodal distribution of RNA expression levels in human skeletal muscle tissue

Mason, C. C. *et al.* *BMC Genomics* 12(), NPG, 2011 [PubMed](#)

Genome-Wide Association Study of Coronary Heart Disease and Its Risk Factors in 8,090 African Americans: The NHLBI CARE Project

Lette, G. *et al.* *Plos Genetics* 7(2), NPG, 2011 [PubMed](#)

Genome-Wide Association Studies of the PR Interval in African Americans

Smith, J. G. *et al.* *Plos Genetics* 7(2), NPG, 2011 [PubMed](#)

Identification and Characterization of an Inborn Error of Metabolism Caused by Dihydrofolate Reductase Deficiency

Banka, S. *et al.* *American Journal of Human Genetics* 88(2), 216-25, 2011 [PubMed](#)

Tumor Growth Inhibition by Olaparib in BRCA2 Germline-Mutated Patient-Derived Ovarian Cancer Tissue Xenografts

Kortmann, U. *et al.* *Clinical Cancer Research* 17(4), 783-91, 2011 [PubMed](#)

Novel Chemosensitive Single-Nucleotide Polymorphism Markers to Targeted
Regimens in Metastatic Colorectal Cancer

Kim, J. C. *et al.* Clinical Cancer Research 17(5), 1200-9, 2011 [PubMed](#)

Ancestry and pharmacogenomics of relapse in acute lymphoblastic leukemia

Yang, J. J. *et al.* Nature Genetics 43(3), NPG, 2011 [PubMed](#)

Acetaminophen-NAPQI Hepatotoxicity: A Cell Line Model System Genome-Wide
Association Study

Moyer, A. M. *et al.* Toxicological Sciences 120(1), 33-41, 2011 [PubMed](#)

Copy number variation and selection during reprogramming to pluripotency

Hussein, S. M. *et al.* Nature 471(7336), NPG, 2011 [PubMed](#)

Pharmacokinetic, pharmacodynamic, and pharmacogenetic determinants of
osteonecrosis in children with acute lymphoblastic leukemia

Kawedia, J. D. *et al.* Blood 117(8), 2340-7, 2011 [PubMed](#)

A genome-wide association study identifies novel loci associated with circulating
IGF-I and IGFBP-3

Kaplan, R. C. *et al.* Human Molecular Genetics 20(6), 1241-51, 2011 [PubMed](#)

Ontogenetic De Novo Copy Number Variations (CNVs) as a Source of Genetic
Individuality: Studies on Two Families with MZD Twins for Schizophrenia

Maiti, S. *et al.* Plos One 6(3), NPG, 2011 [PubMed](#)

Genetic Dosage Compensation in a Family with
Velo-Cardio-Facial/DiGeorge/22q11.2 Deletion Syndrome

Alkalay, A. A. *et al.* American Journal of Medical Genetics Part
a 155A(3), 548-54, 2011 [PubMed](#)

Accuracy of CNV Detection from GWAS Data

Zhang, D. *et al.* Plos One 6(1), NPG, 2011 [PubMed](#)

NEK1 Mutations Cause Short-Rib Polydactyly Syndrome Type Majewski

Thie, C. *et al.* American Journal of Human Genetics 88(1), 106-14, 2011 [PubMed](#)

Replication timing-related and gene body-specific methylation of active human genes

Aran, D. *et al.* Human Molecular Genetics 20(4), 670-80, 2011 [PubMed](#)

Novel common copy number variation for early onset extreme obesity on chromosome 11q11 identified by a genome-wide analysis

Jarick, I. *et al.* Human Molecular Genetics 20(4), 840-52, 2011 [PubMed](#)

Genetic and Metabolic Characterization of Insomnia

Ban, H. J. *et al.* Plos One 6(4), NPG, 2011 [PubMed](#)

Exome sequencing reveals a novel mutation for autosomal recessive non-syndromic mental retardation in the TECR gene on chromosome 19p13

Caliskan, M. *et al.* Human Molecular Genetics 20(7), 1285-9, 2011 [PubMed](#)

Inactivating mutations of acetyltransferase genes in B-cell lymphoma

Pasqualucci, L. *et al.* Nature 471(7337), NPG, 2011 [PubMed](#)

Copy Number Variation in Patients with Disorders of Sex Development Due to 46, XY Gonadal Dysgenesis

White, S. *et al.* Plos One 6(3), NPG, 2011 [PubMed](#)

Genetic Classification of Populations Using Supervised Learning

Bridges, M. *et al.* Plos One 6(5), NPG, 2011 [PubMed](#)

Duplication of MAOA, MAOB, and NDP in a patient with mental retardation and epilepsy

Klitten, L. L. *et al.* European Journal of Human Genetics 19(1), 1-2, 2011 [PubMed](#)

Neural Tumor-Initiating Cells Have Distinct Telomere Maintenance and Can be Safely Targeted for Telomerase Inhibition

Castelo-Branco, P. *et al.* Clinical Cancer Research 17(1), 111-21, 2011 [PubMed](#)

Observation and prediction of recurrent human translocations mediated by NAHR between nonhomologous chromosomes

Ou, Z. *et al.* Genome Research 21(1), 33-46, 2011 [PubMed](#)

Report of Two Brothers With Short Stature, Microcephaly, Mental Retardation, and Retinoschisis-A New Mental Retardation Syndrome?

Phadke, S. R. *et al.* American Journal of Medical Genetics Part a 155A(1), 9-13, 2011 [PubMed](#)

A Novel Homozygous COL11A2 Deletion Causes a C-Terminal Protein Truncation With Incomplete mRNA Decay in a Turkish Patient

Kayserili, H. *et al.* American Journal of Medical Genetics Part a 155A(1), 180-5, 2011 [PubMed](#)

Identification of SOX3 as an XX male sex reversal gene in mice and humans

Sutton, E. *et al.* Journal of Clinical Investigation 121(1), 328-41, 2011 [PubMed](#)

Frequent hSNF5/INI1 Germline Mutations in Patients with Rhabdoid Tumor

Bourdeaut, F. *et al.* Clinical Cancer Research 17(1), 31-8, 2011 [PubMed](#)

Massive Genomic Rearrangement Acquired in a Single Catastrophic Event during Cancer Development

Stephens, P. J. *et al.* Cell 144(1), 27-40, 2011 [PubMed](#)

A SOX9 Duplication and Familial 46,XX Developmental Testicular Disorder

Cox, J. J. *et al.* New England Journal of Medicine 364(1), 91-3, 2011 [PubMed](#)

Aberrant global methylation patterns affect the molecular pathogenesis and prognosis of multiple myeloma

Walker, B. A. *et al.* Blood 117(2), 553-62, 2011 [PubMed](#)

Evolution of human BCR-ABL1 lymphoblastic leukaemia-initiating cells

Notta, F. *et al.* Nature 469(7330), NPG, 2011 [PubMed](#)

Genomic analysis of partial 21q monosomies with variable phenotypes

Roberson, E. D. O. *et al.* European Journal of Human Genetics 19(2), 235-8, 2011 [PubMed](#)

ETV6/RUNX1-positive relapses evolve from an ancestral clone and frequently acquire deletions of genes implicated in glucocorticoid signaling

Kuster, L. *et al.* Blood 117(9), 2658-67, 2011 [PubMed](#)

Assessing Matched Normal and Tumor Pairs in NextGeneration Sequencing Studies

Goh, L. A. *et al.* Plos One 6(3), NPG, 2011 [PubMed](#)

UPD Detection Using Homozygosity Profiling With a SNP Genotyping Microarray

Papenhausen, P. *et al.* American Journal of Medical Genetics Part a 155A(4), 757-68, 2011 [PubMed](#)

Genome-wide association analysis and fine mapping of NT-proBNP level provide novel insight into the role of the MTHFR-CLCN6-NPPA-NPPB gene cluster

Del Greco, M. F. *et al.* Human Molecular Genetics 20(8), 1660-71, 2011 [PubMed](#)

Homozygosity mapping with SNP arrays confirms 3p21 as a recessive locus for gray platelet syndrome and narrows the interval significantly

Fabbro, S. *et al.* Blood 117(12), 3430-4, 2011 [PubMed](#)

A Computational Framework Discovers New Copy Number Variants with Functional Importance

Banerjee, S. *et al.* Plos One 6(3), NPG, 2011 [PubMed](#)

Comparative methylome analysis of benign and malignant peripheral nerve sheath tumors

Feber, A. *et al.* Genome Research 21(4), 515-24, 2011 [PubMed](#)

Sequencing a mouse acute promyelocytic leukemia genome reveals genetic events relevant for disease progression

Wartman, L. D. *et al.* Journal of Clinical Investigation 121(4), 1443-53, 2011 [PubMed](#)

Hunter-gatherer genomic diversity suggests a southern African origin for modern humans

Henn, B. M. *et al.* Proceedings of the National Academy of Sciences of the United States of America 108(13), 5154-62, 2011 [PubMed](#)

Integrative Genomics Identifies Molecular Alterations that Challenge the Linear Model of Melanoma Progression

Rose, A. E. *et al.* Cancer Research 71(7), 2561-71, 2011 [PubMed](#)

Network-Guided Analysis of Genes with Altered Somatic Copy Number and Gene Expression Reveals Pathways Commonly Perturbed in Metastatic Melanoma

Valsesia, A. *et al.* Plos One 6(4), NPG, 2011 [PubMed](#)

Performance assessment of copy number microarray platforms using a spike-in experiment

Halper-Stromberg, E. *et al.* Bioinformatics 27(8), 1052-60, 2011 [PubMed](#)

Reprogramming of Embryonic Human Fibroblasts into Fetal Hematopoietic Progenitors by Fusion with Human Fetal Liver CD34(+) Cells

Sandler, V. M. *et al.* Plos One 6(4), NPG, 2011 [PubMed](#)

IL6-STAT3-HIF Signaling and Therapeutic Response to the Angiogenesis Inhibitor Sunitinib in Ovarian Clear Cell Cancer

Anglesio, M. S. *et al.* Clinical Cancer Research 17(8), 2538-48, 2011 [PubMed](#)

De novo MECP2 duplication in two females with random X-inactivation and moderate mental retardation

Grasshoff, U. *et al.* European Journal of Human Genetics 19(5), 507-12, 2011 [PubMed](#)

The 12q14 microdeletion syndrome: six new cases confirming the role of HMGA2 in growth

Lynch, S. A. *et al.* European Journal of Human Genetics 19(5), 534-9, 2011 [PubMed](#)

Pericentric inversion of chromosome 18 in parents leading to a phenotypically normal child with segmental uniparental disomy 18

Kariminejad, A. *et al.* European Journal of Human Genetics 19(5), 555-60, 2011 [PubMed](#)

Sequence variants at CYP1A1-CYP1A2 and AHR associate with coffee consumption

Sulem, P. *et al.* Human Molecular Genetics 20(10), 2071-7, 2011 [PubMed](#)

Euchromatic 9q13-q21 duplication variants are tandem segmental amplifications of sequence reciprocal to 9q13-q21 deletions

Joseph-George, A. M. *et al.* Journal of Medical Genetics 48(5), 317-22, 2011 [PubMed](#)

SNP-based pathway enrichment analysis for genome-wide association studies

Weng, L. J. *et al.* BMC Bioinformatics 12(), NPG, 2011 [PubMed](#)

